




P260 Rolling Stand and Accessories

For Transporting and Storing the PEDIASCAN® and MAXISCAN® Transilluminators



P260 Rolling Stand

The P260 Rolling Stand is an excellent way to store and transport your PEDIASCAN® or MAXISCAN® transilluminator. The smooth gliding duo-casters allow the unit to be transported from storage or a charging station to bedside with ease.

- 10.5" x 12" (26.67 x 30.48 cm) black powder coated steel tray
- 1" x 48" (2.54 x 121.92 cm) Vel-Cro® hold-down strap
- 2 sets of strap adjustment slots
- 24" (60.96) diameter polished aluminum base
- 39.75" (100.97 cm) overall height.
- Sturdy 1.75" (4.45 cm) diameter chrome plated pole
- Standard non-locking 2" (5.08 cm) duo-casters
- Black powder coated cord wrap

Price: \$595.00/each (domestic pricing)



P260-1 Basket (for P260 Rolling Stand)

The available P260-1 Basket is a great way to store the K913 Battery Charger and assures the unit and charger won't become separated. These baskets are made of steel wire and dipped in vinyl for protection. 12" x 8" x 6" (30.48 x 20.32 x 15.24 cm) in size.

Price: \$135.00/each (domestic pricing)



P260-2 Hospital-Grade Power Strip (for use with P260 Rolling Stand)

The P260-2 Hospital Grade Power Strip is especially designed to fit with the P260 Rolling Stand. Having the power strip available on the rolling stand itself allows the K913 Battery Charger to remain plugged in. Simply plug in the power strip to any U.S. style 110V plug and you will have power to the rolling stand. This power strip is not NEC certified.

Price: \$195.00 (domestic pricing)



The P260 Rolling Stand and its accessories work perfectly with the PEDIASCAN® and MAXISCAN® transilluminators.

CALL 1-800-628-3836 FOR MORE INFORMATION.

Direct: 724•864•9350 - Fax: 724•864•7138 - info@sylvanmed.com - www.sylvanmed.com